







terreno	Salmonella	E. coli	Proteus	Enterobacter sakazakii	Enterobacter cloacae
AS					
BGA					
Cromogeno					
VRBG					
VRBL					
XLD					

terreno	Salmonella	E. coli	Proteus	Enterobacter sakazakii	Enterobacter cloacae
Kliger					
TSI					
Ureasi					
BGB					
Indolo					

Escherichia coli

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


ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

Enterobatteri su BGB

ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

Enterobatteri su BGB

ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

Enterobatteri su BGB

ET: enterotoxinogenic, EP: enteropathogenic, EH: enteronemorragic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

Metodo MPN

Inoculo corrispondente a				MPN/g
1 g	0,1 g	0,01 g		
0	0	0		< 0,3
* 0	1	0		0,3
1	0	0		0,4
* 1	0	1		0,7
1	1	0		0,7
* 1	2	0		1,1
2	0	0		0,9
* 2	0	1		1,4
2	1	0		1,5
* 2	1	1		2
2	2	0		2,1
3	0	0		2,3
* 3	0	1		4
3	1	0		4
3	1	1		7
3	2	0		9
3	2	1		15
* 3	2	2		21
3	3	0		20
3	3	1		50
3	3	2		110
3	3	3		>110

ET: enterotoxinogenic, EP: enteropathogenic, EI: enteroinvasive, DA: diffusely adherent, EA: enteroaggregative,

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ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

BGA

Formula	gm/litre
Proteose peptone	10.0
Yeast extract	3.0
Lactose	10.0
Sucrose	10.0
Sodium chloride	5.0
Phenol red	0.08
Brilliant green	0.0125
Agar	12.0
pH 6.9 ± 0.2 @ 25°C	

Colonial Characteristics

Non-lactose/sucrose-fermenting organisms

[Red-pink-white opaque coloured colonies surrounded by brilliant red zones in the agar - most probably salmonella \(but not Salmonella typhi\).](#)

Proteus and Pseudomonas species

[These may grow as small red colonies.](#)

Lactose/sucrose-fermenting organisms (normally inhibited)

[Yellow to greenish-yellow coloured colonies surrounded by intense yellow-green zones in the agar - Escherichia coli or Klebsiella/Enterobacter group.](#)

ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

Cromogeno per salmonella 1/3

SALMONELLA CHROMOGENIC AGAR BASE

CODE: CM1007

Salmonella Chromogenic Medium is a selective and differential agar base the identification of Salmonella species from other organisms in the family Enterobacteriaceae.

Formula	gm/litre
Special Peptone	10.0
Chromogenic mix	28.0
Agar	12.0
pH 7.2 ± 0.2	

SALMONELLA SELECTIVE SUPPLEMENT

Code: SR0194

Vial contents (each vial is sufficient for 500 ml of medium) per vial per litre		
Cefsulodin	6.0 mg	12.0 mg
Novobiocin	2.5 mg	5.0 mg

ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

Cromogeno per salmonella 2/3

Salmonella Chromogenic Medium is designed to identify Salmonella species based on their utilisation of one chromogenic substrate. Their inability to utilise another chromogenic substrate, that most other members of the family Enterobacteriaceae can utilise, enables rapid and reliable identification of Salmonella species. Traditionally, media used to differentiate Salmonella species from other members of the family Enterobacteriaceae depend upon the ability of Salmonella species to produce hydrogen sulphide coupled with their inability to ferment lactose^{2,3}. These are, however, essentially inadequate methods, with a significant number of the 2000 plus species not exhibiting these characteristics. In recent times chromogenic media have been developed for the rapid and more reliable identification of Salmonella.

Salmonella Chromogenic Agar Base combines two chromogens for the detection of Salmonella sp., 5-Bromo-6-Chloro-3-Indolyl caprylate (Magenta-caprylate) and 5-Bromo-4-Chloro-3-Indolyl b-D galactopyranoside (X-gal). X-gal is a substrate for the enzyme b-D-galactosidase. Hydrolysis of the chromogen, Mag-caprylate, by lactose negative Salmonella species results in magenta colonies.

The medium contains bile salts to inhibit the growth of Gram-positive organisms and the addition of the Salmonella Selective Supplement SR0194 is recommended to increase the selectivity of the medium. This uses novobiocin to inhibit Proteus growth and cefsulodin to inhibit growth of Pseudomonads.

ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

Cromogeno per salmonella 3/3

Species	Colony Colour	Colony Diameter	Colony Morphology
Salmonella spp.	Magenta	1.0 mm	Raised, smooth
Salmonella typhi	Magenta	1.0 mm	Raised, smooth
Salmonella paratyphi	Magenta	1.0 mm	Raised, smooth
Salmonella arizonae	Magenta / blue †	1.5 mm	Raised, smooth
Salmonella gallinarum	Magenta	0.75 mm	Raised, smooth
Salmonella indiana	Blue †	1.0 mm	Raised, smooth
Escherichia coli	Blue	1.0 mm	Raised, smooth
Enterobacter spp.	Blue	1.5 mm	Raised, smooth
Klebsiella spp.	Blue	3.0 mm	Raised, mucoid
Citrobacter spp.	Blue	1.5 mm	Raised, mucoid
Proteus spp.	No growth / straw	0.25 mm	-
Pseudomonas spp.	No growth	-	-
Shigella sonnei	Blue	4.0 mm	Undulate
Shigella dysenteriae	Magenta	1.0 mm	Raised

† Colour of colonies is a presumptive identification as it is dependent on enzyme activity. Some strains of Salmonella arizonae and Salmonella indiana can appear as blue colonies. In addition, some Shigella spp. can appear as magenta colonies. Further confirmatory tests are required.

[Citrobacter](#) – [Escherichia](#) – [Enterobacter](#)
[Salmonella](#) – [Proteus](#) – [Serratia](#) - [Pseudomonas](#)

ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

XLD 1/4

Formula	gm/litre
Yeast extract	3.0
L-Lysine HCl	5.0
Xylose	3.75
Lactose	7.5
Sucrose	7.5
Sodium desoxycholate	1.0
Sodium chloride	5.0
Sodium thiosulphate	6.8
Ferric ammonium citrate	0.8
Phenol red	0.08
Agar	12.5
pH 7.4 ± 0.2 @ 25°C	

ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

XLD 2/4

It relies on xylose fermentation, lysine decarboxylation and production of hydrogen sulphide for the primary differentiation of shigellae and salmonellae from non-pathogenic bacteria. Rapid xylose fermentation is almost universal amongst enteric bacteria, except for members of the *Shigella*, *Providencia* and *Edwardsiella* genera. Xylose is thus included in the medium so that *Shigella* spp. may be identified by a negative reaction. *Salmonella* spp. are differentiated from non-pathogenic xylose fermenters by the incorporation of lysine in the medium. *Salmonellae* exhaust the xylose and decarboxylate the lysine, thus altering the pH to alkaline and mimicking the *Shigella* reaction. However, the presence of *Salmonella* and *Edwardsiella* spp. is differentiated from that of shigellae by a hydrogen sulphide indicator. The high acid level produced by fermentation of lactose and sucrose, prevents lysine-positive coliforms from reverting the pH to an alkaline value, and non-pathogenic hydrogen sulphide producers do not decarboxylate lysine. The acid level also prevents blackening by these micro-organisms until after the 18-24 hour examination for pathogens. Sodium desoxycholate is incorporated as an inhibitor in the medium. The concentration used allows for the inhibition of coliforms without decreasing the ability to support shigellae and salmonellae.

ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

XLD 3/4

Colonial Appearances

Organism	Appearance
<i>Salmonella</i> , <i>Edwardsiella</i>	Red colonies with black centres
<i>Shigella</i> , <i>Providencia</i> , H ₂ S-negative <i>Salmonella</i> (e.g. <i>S. paratyphi A</i>)	Red colonies
<i>Escherichia</i> , <i>Enterobacter</i> , <i>Klebsiella</i> , <i>Citrobacter</i> , <i>Proteus</i> , <i>Serratia</i>	Yellow, opaque colonies

Note

Red colonies may occur with some *Proteus* and *Pseudomonas* species.

[Salmonella](#)

[Escherichia](#), [Enterobacter](#), [Citrobacter](#), [Proteus](#)

[Pseudomonas](#)

ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

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ET: enterotoxinogenic, EP: enteropathogenic, EH: enterohemorrhagic, EA: enteroaggregative, EI: enteroinvasive, DA: diffusely adherent

Corporate (UK) ▾

[Products
by Organism](#)
[Products
by Category](#)
[Product
List](#)
[Safety Data
Sheets](#)
[Technical
Support](#)
[Latest
News](#)
[Exhibitions](#)
[Quality
Certificates](#)
[Literature](#)
[Business
Continuity](#)

Material Safety Data Sheet



Click here to download
Material Safety Data Sheets
for this product.

Organisms

Organisms this product
works with:

■ **Escherichia coli**

Dehydrated Culture Media

[close image gallery](#)

Featured Organism: *Klebsiella*

BRILLIANCE E. COLI/COLIFORM AGAR

Code: CM0956

Brilliance™ E. coli/coliform Agar (formerly Chromogenic E. coli/coliform Agar) is a chromogenic medium to aid differentiation between Escherichia coli and other coliforms in cultures produced from food and environmental samples.

Typical Formula*	gm/litre
Chromogenic mix	20.3
Yeast extract	3.0
Peptone	5.0
Lactose	2.5
Sodium chloride	5.0
Di-sodium hydrogen phosphate	3.5
Potassium di-hydrogen phosphate	1.5
Neutral red	0.03
Agar	15.0
pH 7.0 ± 0.2 @ 25°C	

* Adjusted as required to meet performance standards

Directions

Suspend 55.8g of *Brilliance E.coli/coliform agar* in 1 litre of distilled water. Sterilise by autoclaving at 121°C for 15 minutes. Cool to approximately 50°C. Mix well and pour and into sterile Petri dishes.

Description

Brilliance E.coli/coliform agar is a differential agar used for the presumptive identification of *Escherichia*

coli and coliforms from food and environmental samples. The agar base uses two enzyme substrates to differentiate between *Escherichia coli* and other coliforms. One chromogenic substrate is cleaved by the enzyme glucuronidase, which is specific for *Escherichia coli* and produced by approximately 97% of strains^{1,2,3}. The second chromogenic substrate is cleaved by galactosidase, an enzyme produced by the majority of coliforms. This results in purple *Escherichia coli* colonies, as they are able to cleave both chromogenic substrates, and pink coliform colonies, as they are only able to cleave the galactosidase chromogen.

Technique

Dry the surface of the medium. Prepare the food sample by diluting as appropriate with 0.1% sterile peptone water (CM0009) and homogenise in a stomacher or laboratory blender. Pipette 0.5ml or 1.0ml of the homogenate, as appropriate, on to the plate, and spread over the surface with a glass spreader. Incubate the plates for 18-24 hours at 37°C.

Calculate the total number of coliforms per gram by multiplying purple and pink colonies by the dilution factor. The number of presumptive *Escherichia coli* is obtained by multiplying the number of purple colonies by the dilution factor.

Storage conditions and Shelf life

Store the dehydrated medium at 10-30°C, and use before the expiry date on the label.
Store the prepared medium at 2-8°C.

Appearance

Dehydrated medium; Straw to straw/pink free flowing powder
Prepared medium; Opaque pink coloured gel.

Positive controls:

Escherichia coli ATCC®25922*

Klebsiella pneumoniae ATCC®13883*

Negative control:

Pseudomonas aeruginosa ATCC®27853*

Expected results

Good growth; purple colonies

Good growth; pink colonies

Good growth; straw colonies

* This organism is available as a Culti-Loop®.

Precautions

Refer to the material safety data sheet before handling this product.

Some pathogenic strains, typically *Escherichia coli* O157:H7 are usually glucuronidase negative and therefore will not give purple colonies on this medium.

References

1. Kilian M. and Bulow P. (1976). *Acta. Pathol. Microbiol. Scand.* Sect. B 84: 245-251.
2. Kilian M. and Bulow P. (1979). *Acta. Pathol. Microbiol. Scand.* Sect. B 87: 271-276.
3. Frampton E.W., Restaino L. and Blaszkowski N. (1988). *J. Food Prot.* 51(5): 402-404.

©2001 - 2015 Oxoid Limited, All rights reserved.

Copyright, Disclaimer and Privacy Policy | Conditions of Sale | About Us | Cookies

Thermo Fisher Scientific Inc.

Corporate (UK) ▾

Products
by OrganismProducts
by CategoryProduct
ListSafety Data
SheetsTechnical
SupportLatest
News

Exhibitions

Quality
Certificates

Literature

Business
Continuity

Material Safety Data Sheet

Click here to download
Material Safety Data Sheets
for this product.

Organisms

Organisms this product
works with:

■ **Escherichia coli**

For this Organism

Other products used in the
isolation of [Escherichia
coli](#):

- BR0050
Basic Fuchsin Indicator
- DD0029
Cefpodoxime Combination
Disc Kit
- BR0071
MUG Supplement
- CM0031
Brilliant Green Bile 2%
Broth
- CM0043
MRVP Medium (Clarks and
Lubs Medium)
- CM0451
Lauryl Tryptose Broth
(Lauryl Sulphate Broth)
- CM0505
MacConkey Broth Purple
(US formulation)
- CM0595
Tryptone Bile Agar
- CM0607
Minerals Modified
Glutamate Base
- CM0853
E. C. Broth
- CM0945
TBX Agar
- CM0956
Brilliance E. coli/coliform
Agar
- CM0967
Modified Lauryl Tryptose
Broth with MUG and added
Tryptophan
- CM0968
VIOLET RED BILE
LACTOSE AGAR (ISO)

Dehydrated Culture Media

[close image gallery](#)

Featured Organism: *Escherichia coli*

TRYPTONE BILE X-GLUCURONIDE MEDIUM (TBX)

Code: CM0945

A selective, chromogenic medium for the detection and enumeration of Escherichia coli in food.

Formula	gm/litre
Tryptone	20.0
Bile Salts No. 3	1.5
Agar	15.0
X-glucuronide	0.075
pH 7.2 ± 0.2	

Directions

Suspend 36.6 g of TBX Medium in 1 litre of distilled water. Sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C and pour the medium into sterile Petri dishes.

Description

TBX Medium is based on Tryptone Bile Agar CM0595. Tryptone Bile Agar was originally formulated to improve on earlier methods used to detect *Escherichia coli* in foods^{1,2} in terms of speed, reliability, better recovery from frozen samples and the detection of poor lactose fermenters. TBX Medium builds on these advantages through the addition of a chromogenic agent - X-glucuronide - which detects glucuronidase activity. This is the same enzyme detected by MUG reagent³, and has been shown to be highly specific for *Escherichia coli*⁴. However, approximately 3-4% of *Escherichia coli* are glucuronidase negative, notably *Escherichia coli* O157 strains⁵. Unlike MUG, where the fluorophore leaches out of the cell into the surrounding agar, the released chromophore in TBX Medium is insoluble and accumulates within the cell. This ensures that coloured target colonies are easy to identify. Most *Escherichia coli* strains can be differentiated from other coliforms by the presence of the enzyme

CM0979
E. C. Broth with MUG

CM0980
Lauryl Tryptose Broth with MUG

CM1031
Membrane Lactose Glucuronide Agar (MLGA)

CM1046
Brilliance E. coli/coliform Selective Agar

CM1115
Enterobacteria Enrichment Broth-Mossel (EP/USP/JP/BP)

CM1153
MRS (ISO) AGAR (DE MAN, ROGOSA and SHARPE)

LT0620
Foodborne Pathogens, Monograph Number 5 - E. coli Shigella

Q7085C
Escherichia coli ATCC 8739

DR0300
Dryspot Seroscreen

TD0920
VET-RPLA

TD0960
VTEC-RPLA

QB3648
BAX® REAL-TIME ESCHERICHIA COLI O157:H7 ASSAY

BO0394
MEMBRANE LAURYL SULPHATE BROTH

PO0148
MACCONKEY AGAR WITHOUT SALT

PO0149
MACCONKEY AGAR

PO0727
TBX MEDIUM

PO0745
BRILLIANCE E. COLI / COLIFORM AGAR

PO1142
MACCONKEY AGAR (EP,USP,JP,BP)

PO1226
Brilliance CRE Agar

glucuronidase. The chromogen in TBX Medium is 5-bromo-4-chloro-3-indolyl-beta-D-glucuronide (X-glucuronide), and is targeted by this enzyme. *Escherichia coli* cells are able to absorb this complex intact and intracellular glucuronidase splits the bond between the chromophore and the glucuronide. The released chromophore is coloured and builds up within the cells, causing *Escherichia coli* colonies to be coloured blue/green.

Technique

Dry the surface of the medium in the prepared plates. Dilute the food sample according to the method being followed e.g. 1:10 with Maximum Recovery Diluent CM0733. Homogenise in a stomacher or a laboratory blender.

The following incubation techniques may be used (consult the relevant standard for the complete method):

1. Pipette 0.1 ml of the homogenate on to the plate and spread over the surface with a sterile glass spreader. Incubate the plates for 24 hours at 37°C⁶.
 2. Pipette 0.5 ml of the homogenate on to the plate and spread over the surface with a sterile glass spreader. Incubate the plates for 4 hours at 30°C, then 18-24 hours at 44°C⁷.
 3. Place a cellulose membrane on to the surface of a Minerals Modified Glutamate Medium CM0607 prepared plate. Pipette 1 ml of the homogenate on to the membrane. Incubate for 4 hours at 37°C. Transfer the membrane to a TBX prepared plate and incubate for 18-24 hours at 44°C⁸.
 4. Pipette 1 ml of the homogenate into a sterile Petri dish. Add TBX Medium, cooled to 45°C. Mix well and allow to set. Incubate for 18-24 hours at 44°C. If the presence of stressed cells is suspected pre-incubate the plates for 4 hours at 37°C⁹.
- Multiply the numbers of blue/green colonies by the dilution factor and express the result as the number of *Escherichia coli* per gram of food.

Storage conditions and Shelf life

Store the dehydrated medium at 10-30°C and use before the expiry date on the label.

Store the prepared plates of medium at 2-8°C.

Appearance

Dehydrated Medium: Straw coloured, free-flowing powder.

Prepared medium: Straw coloured gel.

Quality control

Positive controls:

Escherichia coli ATCC® 25922 *

Negative control:

Klebsiella pneumoniae ATCC® 29665

* This organism is available as a Culti-Loop®

Expected results

Good growth; blue/green coloured colonies.

Good growth; straw coloured colonies.

References

1. Gross R.J. and Rowe B. (1985) *J. Hyg. Camb.* 95. 513-550.
2. Anderson J.M. and Baird-Parker A.C. (1975) *J. Appl. Bact.* 39. 111-117.
3. Feng P.C.S. and Hartmann P.A. (1982) *Appl. Environ. Microbiol.* 43. 1320-1329.
4. Hansen W. and Yourassowsky E. (1984) *J. Clin. Microbiol.* 20. 1177-1179.
5. Ratnam S., March S.B., Almed R., Bezanson G.S. and Kasatiya S. (1988) *J. Clin. Microbiol.* 26. 2006-2012.
6. Donovan T.J. *et al* (1998) *Communicable Disease and Public Health* 1 : 188-196.
7. PHLS Standard Methods for Microbiological Examination of Food, Dairy and Water Samples. F20: Direct Enumeration of *Escherichia coli*.
8. ISO 16649-1: 2001. Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli*. Part 1: Colony-count technique at 44°C using membranes and 5-bromo-4-chloro-3-indolyl-beta-D-glucuronide.
9. ISO 16649-2: 2001. Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli*. Part 2: Colony-count technique at 44°C using 5-bromo-4-chloro-3-indolyl-beta-D-glucuronide.

Corporate (UK) ▾

Products by Organism

Products by Category

Product List

Safety Data Sheets

Technical Support

Latest News

Exhibitions

Quality Certificates

Literature

Business Continuity

Material Safety Data Sheet

Click here to download Material Safety Data Sheets for this product.

Optional Products

SR0172
C-T Supplement

Organisms

Organisms this product works with:

Escherichia coli O157

For this Organism

Other products used in the isolation of [Escherichia coli O157](#):

CM0813
Sorbitol MacConkey Agar

CM0981
Sorbitol MacConkey Agar (SMAC) with BCIG

CM0989
Modified Tryptone Soya Broth

CM0990
E. C. Broth (Reduced Bile Salts)

CM1005
Cefixime Rhamnose Sorbitol MacConkey Agar Base

SR0172
C-T Supplement

SR0181
Novobiocin Supplement

SR0191
Cefixime Supplement

DR0120
Dryspot E. coli O157



DR0300
Dryspot Seroscreen

DR0620
Escherichia coli O157 Latex Test

PO0232
SORBITOL MACCONKEY AGAR

PO0702
SORBITOL MACCONKEY AGAR WITH CEFIXIME TELLURITE

Dehydrated Culture Media

[close image gallery](#)

Featured Organism: *Escherichia coli* O157

SORBITOL MacCONKEY AGAR

Code: CM0813

a selective and differential medium for the detection of Escherichia coli O157

Typical Formula*

	gm/litre
Peptone	20.0
Sorbitol	10.0
Bile salts No.3	1.5
Sodium chloride	5.0
Neutral red	0.03
Crystal violet	0.001
Agar	15.0
pH 7.1 ± 0.2 @ 25°C	

* Adjusted as required to meet performance standards

CEFIXIME-TELLURITE SUPPLEMENT

Code: SR0172

a freeze-dried supplement for use with Sorbitol MacConkey Agar for the selective isolation of Escherichia coli O157

Vial contents	per vial	per litre
Potassium tellurite	1.25mg	2.5mg

Directions

Suspend 51.5g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Allow to cool to 50°C. Pour into sterile Petri dishes.

If required, reconstitute one vial of C-T supplement per 500ml medium, according to instruction in the product leaflet. Aseptically add this to cooled medium before pouring into sterile Petri dishes

Description

Escherichia coli O157 is recognised as a cause of haemorrhagic colitis, an illness characterised by bloody diarrhoea and severe abdominal pain, and haemolytic uraemic syndrome (HUS)^{1, 2, 3, 4, 5}, and as such, it a significant human pathogen.

Sorbitol MacConkey Agar is recommended for the isolation of pathogenic *E. coli* O157. The formulation, based on that described by Rappaport and Henig⁶, is identical to MacConkey Agar No.3, except that lactose has been replaced with sorbitol. *E. coli* O157 does not ferment sorbitol and, therefore, produces colourless colonies. In contrast, most *E. coli* strains ferment sorbitol and form pink colonies. The efficiency of Sorbitol MacConkey Agar has been confirmed by March and Ratnam⁷. These workers reported a sensitivity of 100% and a specificity of 85%, and recommended the medium as a simple, inexpensive, rapid and reliable means of screening for *E. coli* O157.

Chapman and co-workers⁸ added cefixime and potassium tellurite to Sorbitol MacConkey Agar to improve the selectivity of the medium. The level of potassium tellurite selects serogroup O157 from other *E. coli* serogroups and inhibits *Providencia* and *Aeromonas* species. Cefixime is inhibitory to *Proteus* spp.

The use of cefixime and tellurite in Sorbitol MacConkey Agar for isolation of *E. coli* O157:H7 is described in the FDA Bacteriological Analytical Manual⁹.

Technique

1. Make up the agar according to the directions and pour into plates. If necessary dry the surface of the agar.
2. Inoculate the plates with a suspension of the test substance (food, faeces, etc) to produce separated colonies.
3. Incubate at 35°C for 24 hours. Doyle and Schoeni¹⁰ reported that 35-37°C is the optimal temperature for growth of *Escherichia coli* O157. At 44-45.5°C growth is poor, even after 48 hours incubation.

Delay in reading plates beyond 24 hours should be avoided because the colour intensity of sorbitol-fermenting colonies fades, reducing the contrast with non-fermenting colonies. Other Gram-negative organisms including *Pseudomonas*, *Proteus* and *Klebsiella* species are able to grow on Sorbitol MacConkey Agar but may generally be differentiated by the appearance of their colonies.

A diagnostic reagent *Escherichia coli* O157 latex test (DR0620) is available so that instant confirmatory tests can be made from suspicious colonies.

Storage conditions and Shelf life

Dehydrated Sorbitol MaConkey Agar must be stored tightly capped in the original container at 10-30°C.

Cefixime-Tellurite Supplement (SR0172) should be stored in the dark at temperatures below 0°C.

Oxoid Sorbitol MacConkey and Cefixime Tellurite Sorbitol MacConkey plates should be stored in the original packaging, at the temperature stated on the pack or product specification, and protected from direct light. When stored as directed, the unopened product will remain stable until the expiry date on the label.

Locally prepared media can be stored for up to 2 weeks when made from CM0813 and SR0172 according to the manufacturer's instructions and stored at 2-8°C, out of direct sunlight and protected from dessication. A longer shelf life may be attainable, but should be validated under the relevant, local manufacturing and storage conditions.

Appearance

Dehydrated medium: Straw/pink coloured, free-flowing powder
Prepared medium: Dark red gel

Quality control

Unsupplemented

Positive control:

Expected results

Escherichia coli O157:H7 Non-toxigenic NCTC12900 * Good growth; 1-2mm straw colonies

Negative control:

Escherichia coli ATCC® 25922 * Good growth; 1-2mm pink colonies

Supplemented with SR0172

Positive control:

Escherichia coli O157:H7 Non-toxigenic NCTC12900 * Good growth; 1-2mm straw colonies

Negative control:

Escherichia coli ATCC® 25922 * No growth/pinpoint-0.25mm pink colonies

* This organism is available as a Culti-Loop®

Precautions

Although the great majority of *Escherichia coli* O157 strains have a typical appearance on Sorbitol MacConkey Agar, some strains are atypical¹¹. Sorbitol MacConkey Agar cannot be used solely to detect VTEC strains of *Escherichia coli* as some non-toxic strains will not ferment sorbitol¹².

Sorbitol MacConkey Agar and Cefixime Tellurite Sorbitol MacConkey Agar are *in vitro* diagnostic use only, by experienced microbiologists. They must not be used beyond the stated expiry date, or if the product shows any sign of deterioration.

Sterilise specimens, equipment and media properly after use.

Limitations

Micro-organisms with atypical enzyme patterns may give anomalous results.

References

1. Centers for Disease Control 1985 - *United States, 1984, Morbid Mortal Weekly Rep.*, 34. 20-21.
2. Karmali M.A., Petric M., Lim C., Fleming P. C., Arbus G. S. and Lior H. (1985) *J. Infect. Dis.* 151. 775-782.
3. Karmali M.A., Steele B.T., Petric M. and Lim C. (1983) *Lancet i.* 619-620.
4. Pai C. H., Gordon R., Sims H. V. and Bryant L. E. (1984) *Ann. Intern. Med.* 101. 738-742.
5. Waters J. R. (1985) *Can. Dis. Weekly Rep.* 11. 123-124.
6. Rappaport F. and Henig E. (1952) *J. Clin. Path.* 5. 361.
7. March S. B. and Ratnam S. (1986) *J. Clin. Microbiol.* 23. 869-872.
8. Zadik P.M., Chapman P.A. and Siddons C.A. (1993) *J. Med. Microbiol.* 39. 155-158.
9. Food and Drug Administration (1995) *Bacteriological Analytical Manual. 8th Edition. AOAC International. Gaithersburg, MD. Chapter 4, 20-23.*
10. Doyle M. P. and Schoeni S. L. (1984) *Appl. and Envir. Microbiol.* 48. 855-856.
11. Karmali M. A. (1988) *Culture* 9. 2.
12. Lior H. and Borczyk A. (1987) *Lancet. i.* 333.

ATCC® is a registered trademark of American Type Culture Collection.